

2 in 1

ULTRASONIC SCALING × POWDER MAINTENANCE



For oral hygiene specialists.

Varios Combi Pro

COMPLETE ORAL CARE

Ask
your NSK
specialist
for a
demo!



Varios Combi Pro BUNDLE PACK



1 Varios Combi Pro Basic Set comprising:

- PROPHY-Handpiece with 60° Nozzle • PROPHY-Hose
- PROPHY-CHAMBER • FLASHpearl (5 pieces x 15g sachets)
- Maintenance kit • VA2LUX-Handpiece • VA2LUX-Hose • 3 Scaling tips
- Sterilisation case • Irrigation bottle (400mL) • Foot control



1 Varios Combi Pro Perio Set comprising:

- PERIO-CHAMBER • PERIO-Hose • PERIO-Handpiece
- Perio-Mate Nozzle tip • Perio-Mate Nozzle tip remover
- Perio-Mate powder



+ 1 EXTRA
Varios2 VA2-LUX-HP Scaler
Optic Handpiece



+ 1 EXTRA
Perio Handpiece



+ 2 EXTRA
Propy Handpieces



+ 1 EXTRA
Perio-Mate Disposable Nozzle
Tip (Pack of 40)



+ 1 EXTRA
Box of FLASHpearl
(4 bottles)



+ 1 EXTRA
Box of Perio-Mate
Powder (2 bottles)

**BUY
BUNDLE
GET FREE
iCart**



+ FREE
iCart Duo for VCP
with Water Bottle Set

MODEL	ORDER CODE	SALE PRICE	SAVE
Varios Combi Pro Bundle offer	VCP BUNDLPKG	\$15,070	\$1,180
Varios Combi Pro iCart	S9090V	FREE	\$3,765

**Total
Savings
\$4,945**

Varios Combi Pro

Varios Combi Pro Basic Set

- PROPHY-Handpiece with 60° Nozzle • PROPHY-Hose
- PROPHY-CHAMBER • FLASHpearl (5 pieces x 15g sachets)
- Maintenance kit • VA2LUX-Handpiece • VA2LUX-Hose
- 3 Scaling tips • Sterilisation case • Irrigation bottle (400mL)
- Foot control

MODEL	ORDER CODE	SALE PRICE	SAVE
Varios Combi Pro Basic Set	Y1004194	\$7,428 <small>RRP \$9,285</small>	\$1,857



Varios Combi Pro Perio Set

Additional to the Basic set, the Perio Set allows safe subgingival treatment for a range of procedures.



Varios Combi Pro Perio Set

- PERIO-CHAMBER • PERIO-Hose • PERIO-Handpiece
- Perio-Mate Nozzle tip • Perio-Mate Nozzle tip remover
- Perio-Mate powder

MODEL	ORDER CODE	SALE PRICE	SAVE
Varios Combi Pro Perio Set	Y1003042	\$2,853 <small>RRP \$3,356</small>	\$503



The powder spray gently cleans subgingival areas



Gentle on soft tissue

Separating powder and water during application prevents periodontal damage. Powder flow is directed in all directions to remove subgingival plaque.

Automatic spray force adjustment

Connecting the subgingival chamber automatically lowers the spray pressure to 20% of the supragingival level to guarantee safe operation during subgingival application.



Safe plastic nozzle

Because the nozzle tip is flexible, it can follow tooth contours and is easy to insert into pockets thereby minimising patient discomfort. The material is semi-transparent to help with visibility. The nozzle tips are sterile and individually packaged so they are ready to use.

Wide range of applications

Applications range from quick biofilm removal to cleaning around dental implants and periodontal pockets. With the nozzle tip removed, clinicians can remove plaque from subgingival pockets of up to 3 mm.

Biofilm removal inside periodontal pockets



Interdental and gingival margin plaque removal



High Speed Handpieces

S-Max and Ti-Max Series



MODEL **Ti-Max Z800L**
ORDER CODE P1110

- Optic
- 23W Power
- Miniature head

RRP	BUY 1	BUY 2	BUY 3
\$2,342	\$1,991		
		\$1,874	\$1,756



MODEL **Ti-Max Z900L**
ORDER CODE P1111

- Optic
- 26W Power
- Standard head

RRP	BUY 1	BUY 2	BUY 3
\$2,342	\$1,991		
		\$1,874	\$1,756



MODEL **S-Max M800L**
ORDER CODE P1253

- Optic
- 23W Power
- Miniature head

RRP	BUY 1	BUY 2	BUY 3
\$1,633	\$1,388		
		\$1,306	\$1,225



MODEL **S-Max M900L**
ORDER CODE P1254

- Optic
- 26W Power
- Standard head

RRP	BUY 1	BUY 2	BUY 3
\$1,633	\$1,388		
		\$1,306	\$1,225



MODEL **Ti-Max X450L**
ORDER CODE P1078

- Optic
- Rear-exhausted 45° head
- 21W Power

RRP	BUY 1	BUY 2	BUY 3
\$1,982	\$1,685		
		\$1,586	\$1,486



MODEL **S-Max Pico**
ORDER CODE P1140

- Optic
- 8W Power
- Ultra Mini Head

RRP	BUY 1	BUY 2	BUY 3
\$1,633	\$1,388		
		\$1,306	\$1,225



Contra Angle Handpieces

S-Max and Ti-Max Series



MODEL **Ti-Max Z95L**
ORDER CODE C1034SW

- Optic
- 1:5 Speed Increasing
- Switch Spray function
- For FG Burs

RRP	BUY 1	BUY 2	BUY 3
\$2,507	\$2,131		
		\$2,006	\$1,880



MODEL **Ti-Max Z45L**
ORDER CODE C1064

- Optic
- 1:4.2 Speed Increasing
- Switch Spray function
- 45° Head
- For FG burs

RRP	BUY 1	BUY 2	BUY 3
\$2,675	\$2,274		
		\$2,140	\$2,006



MODEL **Ti-Max Z85L**
ORDER CODE C1062

- Optic
- 1:5 Speed Increasing
- Single Spray
- For FG Burs

RRP	BUY 1	BUY 2	BUY 3
\$2,507	\$2,131		
		\$2,006	\$1,880



MODEL **Ti-Max Z25L**
ORDER CODE C1038

- Optic
- 1:1 Direct Drive
- Single spray
- For CA burs (ø2.35)

RRP	BUY 1	BUY 2	BUY 3
\$1,763	\$1,499		
		\$1,410	\$1,322



MODEL **Ti-Max Z24L**
ORDER CODE C1126001

- Optic
- 1:1 Direct Drive
- Single spray
- For CA burs (ø2.35)

RRP	BUY 1	BUY 2	BUY 3
\$1,763	\$1,499		
		\$1,410	\$1,322



MODEL **S-Max M95L**
ORDER CODE C1023

- Optic
- 1:5 Increasing
- Quattro Spray
- For FG burs (ø1.6)

RRP	BUY 1	BUY 2	BUY 3
\$2,084	\$1,771		
		\$1,667	\$1,563



MODEL **S-Max M25L**
ORDER CODE C1024

- Optic
- 1:1 Direct Drive
- Single Spray
- For CA burs

RRP	BUY 1	BUY 2	BUY 3
\$1,599	\$1,359		
		\$1,279	\$1,199



MODEL **S-Max M15L**
ORDER CODE C1025

- Optic
- 4:1 Reducing
- Single Spray
- For CA burs

RRP	BUY 1	BUY 2	BUY 3
\$1,947	\$1,655		
		\$1,558	\$1,460

Non-optic options available. Other connections offered for KaVo, Sirona, W&H, Bien Air.

Guaranteed protection. Certified excellence.

The industry leader in infection control for 30 years.

Our proprietary, best-in-class Clean Head System prevents contaminants from entering the handpiece, while our advanced Non-Retraction Valve protects against cross-contamination. NSK sets the standard for safety and protection with the world's most technologically advanced, award-winning handpiece design.



Ti-Max Nano Handpieces



MODEL **nano95LS**
ORDER CODE C1099

- 1:5 Increasing
- Optic
- Quattro spray
- For FG burs (ø1.6)
- Max speed: 200,000 min⁻¹

RRP	BUY 1	BUY 2	BUY 3
\$2,257	\$1,918		
		\$1,806	\$1,693



MODEL **nano25LS**
ORDER CODE C1100

- 1:1 Direct Drive
- Optic
- Single spray
- For FG burs (ø2.35)
- Max speed: 40,000 min⁻¹

RRP	BUY 1	BUY 2	BUY 3
\$1,586	\$1,348		
		\$1,269	\$1,189



MODEL **nano15LS**
ORDER CODE C1101

- 4:1 Reduction
- Optic
- Single spray
- For FG burs (ø2.35)
- Max speed: 10,000 min⁻¹

RRP	BUY 1	BUY 2	BUY 3
\$2,254	\$1,916		
		\$1,803	\$1,690

A totally new installation style only from NSK



Specification comparison

		Our existing model	NLZ • NLZ E
Motor specs	Number of rotations	2,000-40,000 min ⁻¹	100-40,000 min ⁻¹
	Max torque	3.4 Ncm	4.2 Ncm 120% High Torque
Corresponding model to endodontic mode		—	O (NLZ E)

NLZ-E

Electric Micromotor

PORTABLE SYSTEM WITH
ENDO FUNCTION

+ FREE NLZ Endo Handpiece

- Easily convert air driven handpieces to electric with the all new NLZ-E.
- Correctly transmit reciprocating actions for a smooth endodontic treatment.
- New NLZ Endo 6:1 reciprocating handpiece.
- Endo torque calibration function – This function measures the 'no-load' current values and corrects the torque accordingly, resulting in greatly enhance safety in endodontic treatment.
- Contra-Check Function helps detect abnormalities caused by deterioration of the contra-angle handpiece and insufficient maintenance before use.

MODEL	ORDER CODE	RRP	PROMO
NLZ-E 230 Portable	Y1500361	\$7,023	
NLZ Endo Handpiece Non-optic	C1130	\$1,056	FREE



NLZ E

FREE
with
purchase
of NLZ-E
System



MODEL	ORDER CODE	RRP	PROMO
NLZ S230 Portable	Y1500360	\$5,618	
FX57m Handpiece Optic	C1054001	\$601	FREE

FREE
with
purchase
of NLZ
System



MAINTENANCE

iCare C2/C3

Handpiece oiling and maintenance system

- Pre-set solution volume based on handpiece type
- Thorough cleaning and lubrication with simple operation
- Chuck cleaning and lubricating function
- Purging mode removes excess oil
- Integrated 1.2 litre oil tank reduces burden of refill

MODEL	ORDER CODE	SALE PRICE	SAVE
iCare C2	Y1002787	\$4,865	\$1,622
<ul style="list-style-type: none"> • Standard System Complete • 2 rotation ports • 2 non-rotation ports • 230V • Requires 2 adaptors 			
iCare C3	Y1002788	\$5,197	\$1,732
<ul style="list-style-type: none"> • Standard System Complete • 3 rotation ports • 1 non-rotation port • 230V • Requires 3 adaptors 			



Pana Spray offer

Buy 1 box of 6,
get 1 box of 6
Pana Spray Plus FREE



Pana Spray Plus

Cleaning and lubrication

- For high and low speed handpieces, and air motors
- 500 mL • Alcohol-contained lubricant • Pack of 6 cans

MODEL	ORDER CODE	SALE PRICE	SAVE
Pana Spray Plus Pack	Y900630	\$318	\$318
RRP \$636			

**BUY 1 BOX
GET 1 BOX
FREE**

MOBILE DENTISTRY

VIVA ace

Mobile Dentistry System

Purchase the VIVA ace Mobile Dentistry Unit and get a M95L Redband handpiece, a M25L Blueband handpiece and a iso-E-Lux transillumination diagnostic probe FREE!

MODEL	ORDER CODE	RRP
VIVA ace complete set	Y1004228	\$24,890

S-Max M95L Handpiece

- Optic • 1.5 Increasing • Quattro spray • FG Burs

MODEL	ORDER CODE	RRP	PROMO
S-Max M95L	C1023	\$2,084	FREE

S-Max M25L Handpiece

- Optic • 1.1 Direct Drive • Single spray • CA Burs

MODEL	ORDER CODE	RRP	PROMO
S-Max M25L	C1024	\$1,599	FREE

iso E-Lux Probe

- E-Type Optic Motor Connection • Glass Fiber Optics • Stainless Steel Body

MODEL	ORDER CODE	RRP	PROMO
iso E-Lux Probe	Z1265	\$880	FREE



VIVA ace
MOBILE DENTISTRY SYSTEM

I routinely use my VIVA ace in providing dental services for my patients in Aged Care Facilities or who are housebound. It is lightweight, quiet and is literally Plug and Play. I love It. //

Dr Mark Wotherspoon
BDS (UWA), Principal Dentist -
Dentist To Your Door, CSU - Adjunct
Senior Lecturer, Dr Mark's Hygiene -
Inventor

NSK Oceania Ltd www.nskoceania.co.nz

Product Information | Technical Support

Freephone: 0800 44 33 21 (New Zealand) Email: info@nskoceania.com.au

All offers as listed. Valid until 30 September 2022. Not to be used in conjunction with any other NSK promotions. All prices are as quoted in NZD and exclude GST.